
PEACH BACTERIAL SPOT - DETACHED LEAF ASSAY PROTOCOL (*XANTHOMONAS ARBORICOLA*)

Sterilize the following items using autoclave:

- Paper filter
- Distilled water
- Making media (to keep leaves with humidity during experiment days)
- Bacto™ agar (15g/L) and sterilize by autoclaving. After cool down add propiconazole 2.0 ppm (0.4 ml of 5000 ppm propimax to 1 L media)

INOCULUM PREPARATION:

The strains *X.a.88* (used in this test) and *X.a.42* were isolated from peach. Strains were isolated in glycerol 25% under -80°C conditions in an ultra-low freezer.

X.A. STOCK PREPARATION:

1. Prepare liquid media Difco™ LB Broth, Miller (Luria-Bertani) 25 g in 1000 ml (1L) and sterilize by autoclaving.
2. Use 15 ml tube and add around 8 ml LB media and take a loopful of bacterial growth suspended in sterile LB media and leave for 48 hours to growth in shaker agitation with temperature around 30°C.
3. Spin down 400 µl liquid culture (grown in LB) in a sterile 1.5 ml Eppendorf tube
4. Decant supernatant
5. Add 400 µl 25% glycerol (sterilize by autoclaving) to the pellet
6. Mix and store at -80°C

INOCULUM SUSPENSION:

Inoculum suspensions were prepared from 48-hr-old cultures on Difco nutrient agar. A loopful of bacterial growth was suspended in sterile phosphate buffered saline (1X PBS) containing 0.05% Tween 20 and adjusted to 3 different concentrations: Low ($A_{600nm} = 0.1$), Medium ($A_{600nm} = 0.2$) and High ($A_{600nm} = 0.4$). Control treatment = 1X PBS + 0.05% Tween 20.

MAKE 1 LITER 10X PBS (PHOSPHATE BUFFERED SALINE):

1. 80 g NaCl
2. 2 g KCl
3. 14 g Na₂HPO₄
4. 2.4 g KH₂PO₄
5. 0.05% Tween 20 (add 0.5 ml/1L)
6. Adjust pH to 7.4
7. Sterilize by autoclaving

PLANT MATERIAL:

Collect shoots from (20 -25 cm long) placed inside sealable plastic bags on ice to transport to the lab. Expanded leaves (3rd to 5th leaf from tip) were removed, cut in half across leaf center and used the distal part of the leaf, sterilized by soaking for 5 min in 10% bleach, and then rinsed with sterile distilled water. Ten (10) leaves were assayed per treatment.

